



## NURSE CALL AUDIO STATION

NCA		
	Article number:	NCA VORN_basic NCA VORN_aux NCA VORN_sensor
	Size:	65 x 65 x 18 mm
	Function:	Interface between digital signalling inputs and a digitized audio stream and the Ninthway radio network
	Standards:	EN300-220-1 EN300-220-2 EN300-220-3 IEEE 802.15.4 EN54-25 IIC

### Specifications

#### Functional diagram



## NURSE CALL AUDIO STATION

NCA															
<b>Description</b>	<p>The NCA combines the Voice Over Radio Network (VORN) capability with the functionality of the NCB, the multipurpose warning station.</p> <p>It is equipped with a low voltage audio amplifier for magnet speakers and with a high voltage audio amplifier for electret speakers.</p> <p>It can either be used as a (spoken word) sounder, as an announcement device or as an intercom device.</p>														
<b>NCB warning</b>	<p>It provides 4 switches, 2 built in and 2 external.</p> <p>It provides 4 LED indicators, 2 built in and 2 external.</p> <p>A change of state of one of the inputs triggers the transmission of a data frame with information for the receiving application program.</p> <p>The NCB comes with a number of pre settable operation modes called flavours that couple the input signals to the indicators.</p> <p>Flavour:</p> <table border="0"> <tr> <td>1. Offuse</td> <td>simple I/O device</td> </tr> <tr> <td>2. Door</td> <td>door/window contact</td> </tr> <tr> <td>3. Indicator</td> <td>wireless side indicator</td> </tr> <tr> <td>4. Alarm</td> <td>independent button warning station</td> </tr> <tr> <td>5. Nsalarm</td> <td>dependant button warning station</td> </tr> <tr> <td>6. Triple pull</td> <td>pull cord warning station</td> </tr> <tr> <td>7. Syncalarm</td> <td>dependent two button warning station to be used with wireless Indicator</td> </tr> </table>	1. Offuse	simple I/O device	2. Door	door/window contact	3. Indicator	wireless side indicator	4. Alarm	independent button warning station	5. Nsalarm	dependant button warning station	6. Triple pull	pull cord warning station	7. Syncalarm	dependent two button warning station to be used with wireless Indicator
1. Offuse	simple I/O device														
2. Door	door/window contact														
3. Indicator	wireless side indicator														
4. Alarm	independent button warning station														
5. Nsalarm	dependant button warning station														
6. Triple pull	pull cord warning station														
7. Syncalarm	dependent two button warning station to be used with wireless Indicator														
<b>VORN</b>	<p>It is possible to maintain an audio stream over the Ninthway High Secure Network and produce a life audio link. For that purpose the NCA houses a codec that can open an audio path from the NTM transceiver to several audio outputs.</p> <p>In reverse an audio path can be set up that conveys audio signals from the built-in microphone to the NTM transceiver.</p> <p>In this way a half duplex audio communication over the Ninthway Radio Network can be set up.</p> <p>The codec samples the audio signal at 16 kHz and compresses it using an ADPCM protocol. The digitized audio signal is transmitted @ 500kpb on the BBN frequency band.</p> <p>Full duplex operation requires the use of two audio bands. It requires the use of a double NCA at both sides and repeater stations need to be equipped with a second BBN transceiver that operates on the second audio frequency.</p> <p>Alternatively the codec has space to store pre-recorded messages that can be</p>														

## NURSE CALL AUDIO STATION

NCA	
	played with a simple command to the NCA.
<b>Warning station connections</b>	<p>Built in pull cord switch.</p> <p>Built in reed switch.</p> <p>2 sets of 2 terminals for 2 external buttons.</p> <p>2 sets of 2 terminals for external indicators.</p> <p>4 UART pins, gnd, 5V, rxd, txd.</p>
<b>Audio connections</b>	<p>Electret microphone with built-in AGC.</p> <p>Output pinheads:</p> <ul style="list-style-type: none"> <li>• 1V pp voltage line output.</li> <li>• Current source line output.</li> <li>• 3W@ 40hm loudspeaker connection (Class D).</li> <li>• 4 pinhead I2S digital audio I/O.</li> <li>• 20 V electret speaker output.</li> <li>• Digital I2S audio output bus.</li> </ul>
<b>Options</b>	<p>Temperature sensor MCP9800. -10 to +85 °C.</p> <p>3D Accelerometer ADXL345.</p> <p>Li-ion auxiliary power supply with built in charge controller.</p>
<b>Jumpers</b>	<p>JP8: input connection for line input socket</p> <p>JP12: headphone connection</p> <p>J15: charge select mode 1-2: current; 2-3 voltage controlled</p>
<b>Parameters</b>	See application note 7; Manual for the nurse call audio station
<b>Indicators</b>	<p>LD1 Green LED</p> <p>LD2 Red LED</p> <p>LD3 external</p> <p>LD4 external</p>
<b>Power supply and current consumption</b>	<p>Supply voltage: 2 Terminals 2.5 – 5 V, 0.25 - 1.5 A. Reverse polarity protected. Feeds two supplies.</p> <p>Primary supply: for NTM, NCB and codec: @ 250 mA.</p> <p>Secondary supply: for built in audio amplifier using same power supply input 5V @1.5A. Can be shut down to save power.</p> <p>Option: Auxiliary power with 500 mA charger, Li-ion A-cell (1200 mAh) when secondary supply present.</p>

## NURSE CALL AUDIO STATION

NCA	
<b>Radio parameters</b>	<p>Audio data is transmitted real time over the network in 80 frames per second using the CSMA-CA protocol. This exceeds the duty cycle regulations for bands like 868 MHz. Therefore VORN operations takes place on the 863 – 865 MHz band.</p> <p>This requires gateway and repeater stations to have their BBN transceiver set to the 863 MHz band and their data rate to 500 kbps.</p> <p>The settings for the sensor and actor network band stays @ 868.3 MHz and 100 kbps data rate. Only during VORN operation will the transceiver on the NCA set itself to the BBN band and data rate.</p>
<b>Audio parameters</b>	<p>The NCA houses four audio outputs.</p>
Line output	<p>The head connection (j10) 1-2-3 provides two single ended line/headphone outputs.</p> <p>1-2 current output: max 6 mA            2-3 voltage output: max 1Vpp            Minimum load impedance: 5 KOhm            DC_bias: 1.2V</p>
Audio amplifier for coil speaker	<p>Mag spkr connection provides BTL coil speaker output.</p> <p>Output power @ 4 Ohm 2.30 W            Output power @ 8 Ohm 1,36 W</p> <p>Harmonic distortion @ 4 Ohm 0.03%            Harmonic distortion @ 8 Ohm 0.02%</p> <p>Fixed amplifier gain: 18 dB            SNR 105 dB            Bandwidth: 0 – 50 KHz.</p>
HV output	<p>HV output, class G ceramic speaker driver.</p> <p>Max output voltage: 20 Vpp            Clip-and-plop level: -67 dBV            Continuous output power: 2.4 W            Maximum capacitive load: 1µF            Fixed amplifier gain: 32 dB            Bandwidth: 0 – 10KHz            SNR 105 dB</p>
<b>Codec audio output</b>	<p>Standard volume attenuation: 12 dB            Minimum attenuation: 0 dB            Maximum attenuation: 63.75 dB            AGC default max gain: 28 dB</p>

## NURSE CALL AUDIO STATION

NCA	
	AGC maximum gain: 47 dB
<b>NTM mode</b>	The NCA-mode is activated using function 7.
<b>Mounting instructions</b>	For the best performance of the radio transceiver mount the PCB with the antenna upright when possible.
<b>Additional information</b>	Datasheet NTM_3 Application note 1: Programming the NTM Application note 2: Ninthway high secure radio network Application note 7: Manual for the nurse call audio station